

Novozymes' capabilities in biocatalysis for pharmaceutical industry

诺维信在制药行业的生物催化解决方案

Novozymes Presentation

November 2015

❑ NOVOZYMES AND BIOCATALYSIS 诺维信和生物催化

❑ CASE STUDIES 实例分析

Pharmaceuticals & fine chemicals are one of the most polluting industries

制药&精细化工是最污染的行业之一

Environmental impact of some of the chemical industries 化工业对环境的若干影响

Industries 行业	Production Scale 生产规模 (Tons per Year) (吨/年)	E-factor E-因子
Oil Refining 炼油业	$10^6 - 10^8$	< 0.1
Bulk Chemicals 大宗化学品业	$10^4 - 10^6$	$< 1 - 5$
Fine Chemicals Industry 精细化工	$10^2 - 10^4$	$5 - > 50$
Pharmaceutical Industry 制药业	$10 - 10^3$	$25 - > 100$

A sustainable and greener footprint in chemical production is increasingly becoming a global need...一条可持续的环保的化工生产道路正日益成为全球所需

There are existing and emerging solutions to counter pollution 这是目前新起的防止污染的解决方案


Biotransformation is one of the key sustainable answer
to growing concerns of chemical industry

生物转化是化工业日益受关注的关键原因



Fermentation

发酵



Biocatalysis (Isolated Enzymes)

生物催化 (隔离酶)



Whole Cells

全细胞

There are several benefits of using enzymes over chemical reaction to produce APIs

生产APIs的化学反应中使用酶制剂有若干好处

General Advantages using Enzymes in API Manufacturing 在APIs生产过程中使用酶制剂的综合优势



- Improved Productivity 提高生产力
 - Higher Yields 更好的产量
 - Shortened synthesis route 缩短合成路线
- Greater Savings Potential 更加节约成本
 - Replaces costly chiral resolving agents 替代价格昂贵的手性拆分剂
 - Significant reduction in waste streams 有效地减少污水排放
- Greater Selectivity 更高的选择性
 - High level of Stereo-, regio-, and chemo-selectivity 高度的立体，区域和化学选择性
- Milder reaction conditions 更温和的反应条件
 - Mostly ambient reaction condition 主要是常温反应条件
- Simplified work-streams 简化工作流程
 - Simplified processing and purification 简化加工和纯化
 - Fewer byproducts, reducing impurities 更少的副产品，减少杂质

Enzymes are already used in producing some of very important fine chemicals

酶制剂已经用于生产一些重要的精细化学品

Some examples of biocatalysis application in producing fine chemicals & API 生物催化在精细化学品&API生产中的应用实例

Compound 化合物	Enzyme 酶制剂	Manufacturer 制造商
Chiral amines 手性的胺类	<i>C. antarctica lipase B</i>	BASF
(R)-Mandelic acid 扁桃酸	<i>(S)- or (R)-nitrilases</i>	BASF
L- or D-Amino acids 氨基酸	<i>amidase</i>	DSM
Chiral alcohol 手性的乙醇	<i>lipase*</i>	DSM
2-halopropionic acids 卤丙酸	<i>porcine lipase</i>	DSM
(R)-Glycidyl 缩水甘油基	<i>ester lipase</i>	DSM
β-Phenylalanines 苯丙氨酸	<i>lipase</i>	Dowpharma
Carbocyclic nucleosides intermediate 碳环型核苷中间体	<i>γ-lactamase</i>	Dowpharma
Protein farnesyl transferase inhibitor 蛋白质西基转氨酶抑制剂	<i>lipase LIP-300</i>	Schering
Pregabalin 普瑞巴林	<i>Lipozyme TL 100 L(Novozyme enzyme)</i>	Pfizer
Carbocyclic nucleosides intermediate 碳环型核苷中间体	<i>Savinase (Novozyme enzyme)</i>	GSK/other players

Most of these reaction are run at **industrial (multi Ton) scale**大部分反应为工业规模（万吨级）

Based on the reaction types, enzymes are grouped into several classes

根据反应类型，酶制剂可分为若干类

Hydrolases are the most used enzyme class in biotransformation's followed by oxidoreductases 水解酶是生物转换中最常用的酶，其次是氧化还原酶

Some key examples of enzyme classes and their applications 酶的种类和应用实例

Selected Reactions 参考反应	Enzyme Class 酶种类	EC No.	Key Applications 核心应用
- C=O and C=C reduction 降低 Reductive amination of C=O 还原胺化反应 Oxidation of C-H, C=C, C-N and C-O 氧化 - Cofactor reduction/oxidation 辅因子降低/氧化	Oxidoreductase 氧化还原酶	1	- Atorvastatin Intermediate阿托法他汀中间体 - Montelukast孟鲁司特 - Sitagliptin西格列汀
Transfer of functional groups such as amino, Acyl, phosphoryl, methyl, glycosyl, nitro & sulphur groups 功能基团的转化，如氨基、酰基、磷酸基、甲基、糖基、硝基&硫磺基团	Transferase 转氨酶	2	Cyclodextrin from Starch 淀粉提取环糊精
Hydrolysis of esters, amides, lactones, lactams, epoxides, nitriles & reverse reactions 酯类，酰胺类，内酯类，内酰胺类，环氧化合物、腈类的水解反应&逆反应	Hydrolase 水解酶	3	- Diltiazem Intermediate 地尔硫卓中间体 - Resolution of rac. nylethylamine 外消旋体拆分
Addition of small molecules to double bonds such as C=C, C=N and C=O 双键上添加小分子	Lyase 裂解酶 (synthase)		Apple acid intermediate他汀
Isomerisations such as racemizations,	Isomerase		HECS from glucose

Hydrolase 水解酶

Novozymes has
 a comprehensive
 Lipase/Protease/
 Amylase
 portfolio诺维信拥有广泛的脂肪酶/蛋白酶/淀粉酶

Novozymes has one of the most comprehensive biocatalysis enzyme portfolio

诺维信拥有最广泛的生物催化酶组合




- Lipases 脂肪酶
- Proteases 蛋白酶
- Esterases/cutinases 酯酶
- Oxido-reductases; 氧化还原酶
 - laccases, 漆酶
 - peroxidases, 氧化物酶
 - Chloroperoxidase 氯化物过氧化酶
- Amylolytic enzymes 淀粉酶

Novozym 435- is the industry gold standard in lipases

NOVOZYM 435-是脂肪酶的行业黄金标准

Novozym 435 is CALBL immobilized on a polymer support

Novozym 435 是将CALBL固定于聚合物载体上



Key Advantages of Novozym 435 主要优势

- Mild and selective on multi-functional substrates 温和和选择性的多功能底物
- Active both in bulk liquid substances and organic co-solvents 在散装液体底物和有机共溶剂均保持活跃
- Functions across wide temperature range (20-110°C) 作用于宽温域
- Suitable for both stirred batch tank / Continuous fixed bed reactors 适用于批次搅拌槽和连续固定床反应器
- Recycled without activity loss for 5-10 or even more times, depending on reaction conditions 根据反应条件，可回收利用5-10次，甚至更多次不失活力
- Large-scale industrial production 大规模工业化生产
- Nil residue in the final products 最终产品无残留
- Works in anhydrous conditions & so can be used for moisture sensitive substrates 脱水条件下作用&可用于湿度敏感性底物
- Can be used with metal based racemizing reagents 可与金属外消旋试剂一起使用

Key applications of Novozym 435 主要应用

- Dynamic kinetic resolution of **amine** coupled with metal racemising agent 动态动力学拆分金属外消旋溶剂配合的胺类
- Dynamic Kinetic resolution of **alcohols** coupled with ruthenium catalyst 动态动力学拆分钌催化剂配合的醇类
- Kinetic resolution of **Darunavir Intermediate** 动力学拆分地瑞纳韦中间体

Novozymes Immobilized Lipases.....

诺维信固定化脂肪酶

Novozym®435

CALB enzyme adsorbed on
Macroporous Acrylic Resin

脂肪酶B吸附于丙烯酸树脂

Lipozyme®TL IM

TL enzyme adsorbed on
Silica Gel

TL酶吸附于硅胶

Novozym 40086 *

RM Enzyme adsorbed on
Macroporous Anionic Resin

RM酶吸附于阴离子型树脂

* Previously it was "Lipozyme ® RM IM"

We have the industry leading robust lipase portfolio

我们拥有行业标杆的脂肪酶组合

CalA type of enzymes

- *candida antarctica* lipase A
南极假丝酵母脂肪酶A
- Allows branching on C-alpha and bulky side chains on the alcohol part
能够在C-alpha 和醇基庞大链一侧分支

Key Applications 主要应用

- Kinetic Resolution of 2-phenylbut-3-yn-2-ol
动力学拆分2-苯基丁-3炔-2-醇
C10H10O
- Resolution of 1-methyl-1,2,3,4-tetrahydronaphthalen-1-ol
拆分 1-甲基-1,2,3,4-四氢萘-1-醇
C11H14O
- Kinetic resolution of (-)-**Paroxetine Intermediate**
动力学拆分帕罗西汀中间体

CalB type of enzymes

- *candida antarctica* lipase B
南极假丝酵母脂肪酶B
- Allows larger groups on acid part and high specificity on alcohol part 允许较大的基团结合酸根并且对醇基有特异结合

Key Applications 主要应用

- Kinetic resolution of **Lotrafiban key intermediate** 动力学拆分洛曲非班重要中间体
- Kinetic Resolution of (S)-(+)-**Citalopram Intermediate** 动力学拆分西酞普兰中间体
- Kinetic Resolution of **Rivastigmine Intermediate** 动力学拆分卡巴拉汀中间体
- Kinetic resolution of **(R,R)-Formoterol intermediate** 动力学拆分福莫特罗中间体
- Kinetic resolution of **Ketoprofen&Flurbiprofen**

TLL type of enzymes

- *Thermomyces lanuginosus* lipase
细毛嗜热霉脂肪酶
- Allows larger groups on both alcohol and acid parts 允许较大基团结合在醇基和酸根上

Key Applications 主要应用

- Kinetic Resolution by desymmetrization of CNDE-key intermediate for **Pregabalin**
通过CNDE重要中间体去对称化动力学拆分普瑞巴林

Some of the other important enzymes in Novozymes portfolio

are the proteases- Alcalase & Savinase

诺维信其他重要的蛋白酶—Alcalase and savinase



Alcalase碱性蛋白酶

- *subtilisin* based alkaline proteases
枯草杆菌碱性蛋白酶
- Available as liquid formulations
Alcalase 2.4LFG & Alcalase 2.5L
有液体剂型Alcalase 2.4 L FG 和 Alcalase 2.5L

Key applications 主要应用

- Dynamic kinetic resolution of protected Amino acids eg. DKR of ***Benzylated Tyrosine***
动力学拆分保护氨基酸 如苄化酪氨酸

Savinase蛋白酶

- *subtilisin* based alkaline proteases
枯草杆菌碱性蛋白酶
- Available as Savinase 12T solid granulate & liquid formulation
Savinase 16L
有固体颗粒Savinase 12T & 液态剂型Savinase 16L

Key applications

- Kinetic resolution of vince lactam-key intermediate for ***Carbovir*** & ***Abacavir***
动力学拆分2-氮杂双环[2.2.1]庚-5-烯-3-酮(文斯酯)—
—卡巴韦和阿巴卡韦的重要中间体

- ❑ NOVOZYMES AND BIOCATALYSIS 诺维信和生物催化
- ❑ CASE STUDIES 实例分析


Pregabalin Biocatalytic Route

普瑞巴林的生物催化路径

Improvements due to replacing chemical process with chemoenzymatic reaction


化学酶法反应替代化学过程所带来的改善

Pregabalin Chemical synthesis 普瑞巴林的化学合成




在外消旋的普瑞巴林阶段是无效拆分

Pregabalin Biocatalytic Route 普瑞巴林的生物催化路径



Flow Chart of Pregabalin Process.....

普瑞巴林流程图



Key Advantages of Biocatalytic Route 生物催化路径的主要优势

- Low protein Loading (0.8%) 低酶蛋白用量
- Resolution at first stage-wrong isomer can be recycled 在第一阶段拆分—错误的同分异构体可回收利用
- High throughput 高流通量
- All reactions conducted in water 所有反应在水中
- E-factor improved from **86** to **17** E-因子从86减少到17

- Starting material reduction by 800 Mt (cumulative) 初始物料减少800 Mt (累计的)
- Drastic reduction in solvent usage 显著减少溶剂的使用
- Mandelic acid usage-500Mt eliminated (Cumulative) 扁桃酸的使用减少 500 Mt (累计的)
- Energy usage reduced by 83% 能源消耗降低 83%

Pregabalin Chemoenzymatic route used >5x less inputs than Chemical Route

化学酶法合成普瑞巴林比化学法少了多于5倍的投入量

Raw Material Inputs for 1000kg API Output

1000公斤API的原料投入量

Inputs 投入	Chemical Route 化学路径 Kg	Chemoenzymatic Route 化学酶法路径 Kg	Reduction 减少率
CNDE	6212	4798	23%
Enzyme	0	574	
S-Mandelic acid	1135	0	85%
Raney Nickel	531	79.5	
Solvents	50042	6230	88%
Total	57920	11681.5	

Trandolapril Biocatalytic Route



群多普利的生物催化路径

Novozym 435 gives optimum performance in stereo selective synthesis of Trandolapril key intermediate

Novozym 435可以高效催化群多普利重要中间体的立体选择性合成

Trandolapril Key Intermediate Synthesis 群多普利重要中间体的合成

Trandolapril Biocatalytic Route 群多普利生物催化路径



Key Advantages of Biocatalytic Route: 生物催化路径的主要优势:

- Enzyme loading is low 低酶加量
- Substrate concentration can be as high as 100g/L 底物浓度高达100g/L
- Yield is high ~80-90% & ee> 98% 得率高达80-90% & ee>98%
- Enzyme can be recycled for approximately 5 times with out fresh loading 酶大概可以回用5次
- Reaction conditions i.e. pH & temperature is not very harsh 反应条件温和
- Reaction can be converted from batch to continuous process 反应可以由批式转为连续过程


Vince Lactam-Carbovir & Abacavir Key Intermediate Biocatalytic Route

VINCE LACTAM-卡巴韦及阿巴卡韦重要中间体生物催化路径


Savinase 12T is used on a commercial scale to produce Vince Lactam- Key Raw Material of Carbovir & Abacavir

Savinase 12T用于商业规模生产Vince Lactam-卡巴韦及阿巴卡韦的重要原料

Carbovir & Abacavir Synthesis 卡巴韦及阿巴卡韦的合成路径



Vince Lactam Biocatalytic Route Vince Lactam生物催化路径



Key Advantages of Biocatalytic Route:生物催化路径的主要优势:

- Substrate concentration can be as high as 100g/L 底物浓度高达100g/L
- Conversion is ~50% & ee of 1R,4S Vince lactam >99% 转化率为50% & ee of 1R, 4S Vince lactam >99%
- Reaction can be carried out in a mixture of solvent & phosphate buffer 反应在有机溶剂和磷酸缓冲液中进行
- Reaction is carried out under mild & ambient conditions 反应条件温和



Posaconazole Key Intermediate Biocatalytic Route

泊沙康唑重要中间体生物催化路径

Novozym 435 gives optimum performance in stereoselective synthesis of Posaconazole key intermediate
Novozym 435可以高效催化泊沙康唑重要中间体的立体选择性合成

Posaconazole Synthesis 泊沙康唑合成路径

Posaconazole Intermediate Biocatalytic Route 泊沙康唑中间体生物催化路径



Key Advantages of Biocatalytic Route:生物催化路径的主要优势:

- Substrate concentration can be as high as 100g/L 底物浓度高达100g/L
- Conversion is 70-80% 转化率为70-80%
- Reaction can be carried out in a non polar solvent like toluene 反应可以在非极性溶剂如甲苯中进行
- Immobilized enzyme can be recycled 固定化酶可被回收再利用



Darunavir Key Intermediate Biocatalytic Route

地瑞那韦重要中间体生物催化路径

Lipozyme CALBL is commercially used for stereoselective synthesis of key bicyclic intermediate for Darunavir
Lipozyme CALBL用于商业规模立体选择性合成地瑞那韦重要中间体

Darunavir Synthesis 地瑞那韦合成路径

Darunavir Intermediate Biocatalytic Route 地瑞那韦中间体生物催化路径



Key Advantages of Biocatalytic Route: 生物催化路径的主要优势:

- Substrate concentration can be as high as 100g/L 底物浓度高达100g/L
- Reaction can be carried out using lipozyme CALBL or Novozym 435 可以使用Lipozyme CALBL 或者Novozym 435
- Reaction can be carried out using water as a solvent 在水中进行
- Reaction is carried out under mild & ambient conditions 反应条件温和


Latanoprost & Travoprost Biocatalytic Route

拉坦前列素和曲伏前列素生物催化路径

Novozym 435 is used for the key regioselective esterification step


Novozym 435用于关键区域选择性酯化步骤

Latanoprost & Travoprost Synthesis 拉坦前列素和曲伏前列素合成路径



Latanoprost & Travoprost Biocatalytic Route

拉坦前列素和曲伏前列素生物催化路径



Key Advantages of Biocatalytic Route: 生物催化路径的主要优势:


- Biocatalytic reaction is carried out at high substrate concentrations. 高底物浓度
- Reaction yields are >90% & highly regioselective. 得率>90%, 高区域选择性
- Reaction can be carried out using appropriate solvent or without any solvent 使用适合的溶剂或是不用溶剂
- Reaction is carried out under mild & ambient conditions 反应条件温和

Moxifloxacin Key Intermediate Biocatalytic Route


莫西沙星重要中间体生物催化路径

Key step for moxifloxacin key intermediate manufacturing is enzymatic
酶可以催化莫西沙星重要中间体生产的关键步骤

Moxifloxacin Side Chain Synthesis 莫西沙星侧链合成路径




Moxifloxacin Synthesis 莫西沙星合成路径



Moxifloxacin Key Intermediate Biocatalytic Route

莫西沙星重要中间体生物催化路径

Lipozyme CALBL or Novozym 435 can be used for desymmetrisation of racemic substrate to the required 2S,3R isomer
Lipozyme CALBL或者Novozym 435用于去对称化，即外消旋底物成2S, 3R异构体




Key Advantages of Biocatalytic Route:生物催化路径的主要优势:

- Biocatalytic reaction is carried out at high substrate concentrations ~100g/L 底物浓度高达100g/L
- Reaction yields are >90% & highly regioselective. 得率>90%,高区域选择性
- Reaction can be carried out using appropriate solvent or without any solvent in water with appropriate buffer 可以在适合的溶剂中进行，也可以在缓冲液中进行
- Reaction is carried out under mild & ambient conditions 反应条件温和


Repaglinide Key Intermediate Biocatalytic Route 瑞格列奈重要中间体生物催化路径

Key step for repaglinide key intermediate manufacturing is enzymatic
酶可以催化瑞格列奈重要中间体生产的关键步骤

Repaglinide Side Chain Synthesis 瑞格列奈侧链合成路径




Repaglinide Synthesis 瑞格列奈合成路径



Repaglinide Key Intermediate Biocatalytic Route

瑞格列奈重要中间体生物催化路径

Novozym 435 is used for regioselective hydrolysis of diester to the required intermediate
Novozym 435用于区域选择性水解二元酯产生必需中间体



Key Advantages of Biocatalytic Route: 生物催化路径的主要优势:


- Biocatalytic reaction is carried out at high substrate concentrations ~200g/L 底物浓度高达200g/L
- Reaction yields are ~100% & highly regioselective. 得率100%, 高度区域选择性
- Reaction can be carried out using appropriate solvent or without any solvent in water with appropriate buffer 反应在适当溶剂中进行, 或是在缓冲液系统中
- Reaction is carried out under mild conditions 反应条件温和

L-Carnitine Key Intermediate Biocatalytic Route


左旋肉碱重要中间体生物催化路径

Key step for L-Carnitine key intermediate manufacturing is enzymatic
酶可以催化左旋肉碱重要中间体生产的关键步骤

L-Carnitine Key Intermediate Synthesis 左旋肉碱重要中间体合成路径




L-Carnitine & Acetyl Carnitine Synthesis
左旋肉碱和乙酰肉碱重要中间体合成路径



L-Carnitine Key Intermediate Biocatalytic Route

左旋肉碱重要中间体生物催化路径

Novozym 435 can be used for kinetic resolution of racemic ethyl-4-Chloro-3-hydroxybutyrate to ethyl-(R)-4-Chloro-3-hydroxybutyrate Novozym 435 用于动力学拆分以消旋4-氯-3-羟基丁酸乙酯成为(R) 4-氯-3-羟基丁酸乙酯





Key Advantages of Biocatalytic Route:生物催化路径的主要优势:

- Biocatalytic reaction is carried out at high substrate concentrations ~100g/L 底物浓度高达100g/L
- Reaction yields are ~70% of theoretical yield & highly stereoselective ~99.9% ee 得率70%，高度立体选择性
- Reaction can be carried out using appropriate solvent or without any solvent in water with appropriate buffer 反应在适当溶剂中进行，或是在缓冲液系统中
- Reaction is carried out under mild & ambient conditions 反应条件温和

Case story: Biocatalysis for cosmetic

Making chemicals more efficiently

Emollient ester




Esterification for production of cosmetic fatty acid esters.
Chem. Soc. Rev., 2013, **42**, 6475--6490


lipases and their hydrolytic, esterifying and acylating activities show enormous potential for implementation in the production of cosmetic ingredients

Lipase Catalysed Reactions


Hydrolysis




Esterification




Interesterification




Transesterification (e.g. methanolysis)




Epoxidation




Polycondensation



Polycondensation



Ring opening





let's discuss how best we could work together

让我们共同探讨如何更好地合作

CONTACT

联系方式

NAME : ZHE HUANG AUSTIN

TITLE : REGIONAL BUSINESS MANAGER

EMAIL : ZHU@NOVOZYMES.COM

NAME : JUNXIN DUAN

TITLE : TECHNICAL SERVICE MANAGER

EMAIL : JDUA@NOVOZYMES.COM